

K'Ball Universal End Pin Protector

K'Ball® is indicated in the treatment of traumas of the upper and lower extremities. It is used in conjunction with the protruding wire to prevent its inward migration.

K'Ball® is a non-reusable device that fits all wires ranging a \varnothing from 0.8 to 2.0 mm. It protects patients and hospital staff from sharp K-wire ends, and shields the surgeon from cuts and pricks when multiple wires are handled in a small area.

K'Ball® is designed to withstand repeated steam sterilizations.

K'Ball® does not require a colour code nor the decried 90° bending procedure of the protruding end.

K'Ball® can be loosen and tighten at will to ease the assembly and removal procedures.

K'Ball® is designed on the principle of simplicity and intuition. In this manner either left or right lateral hole can be used to slip the device onto the protruding end, whatever its diameter.

K'Ball® is an efficient healthcare solution with a minimum inventory requirement that also enables to reduce the volume of dressing.



Distributed By:
Lavender Medical Limited
Unit 28, Wedgwood Way
Steveange, SG14QT

Tel: 0845 676 9733
Email: orders@lavendermedical.com
Web: www.lavendermedical.com