

FOOTMOTION

PERCUTANEOUS SURGERY




NEWCLIP-TECHNICS


► **Application** : FOOTMOTION percutaneous surgery offers a range of instruments dedicated to Forefoot Surgery.

PERCUTANEOUS REAMERS

Ø 2.0 mm/Ø 2.2 mm SHANNON REAMER (straight or helical flute)



Ø 3.0 mm SHANNON REAMER (straight or helical flute)



Ø 4.0 mm WEDGE REAMER



WEDGE CONICAL REAMER



APPLICATIONS

Percutaneous reamers are designed for minimally invasive surgery of the forefoot*:

► **Cylindrical reamers**
- Ø4.0 mm wide and Ø3.0 mm long Shannon reamer: exostosectomy, arthrodesis and shortening osteotomy.

- Ø2.0 mm short and long Shannon reamer: lateral rays osteotomy, distal or proximal osteotomy of the first ray and osteotomy of the first phalanx.

► **Conical reamers**
- Wedge: Distal monocortical osteotomy of the first ray and osteotomy of the first phalanx.

REFERENCES

PERCUTANEOUS REAMERS**

Ref	Description
ANC197	Ø4 mm large Shannon reamer - Helical flute L 16 mm
ANC198	Ø2.2 mm short Shannon reamer - Helical flute L 12 mm
ANC199	Ø4 mm wedge reamer - Flute L 13 mm
ANC203	Ø3 mm wedge reamer - Flute L 13 mm
ANC476	Ø2.0 mm long Shannon reamer - Helical flute L 20 mm
ANC537	Ø2.0 mm short Shannon Isham reamer - Straight flute L 13 mm

PERCUTANEOUS REAMERS**

Ref	Description
ANC538	Ø2.2 mm long Shannon Isham reamer - Straight flute L 22 mm
ANC599	Ø2.0 mm ultra short Shannon reamer - Helical flute L 8 mm
ANC661	Ø3.0 mm long Shannon reamer - Helical flute L 21 mm
ANC842	Ø3.0 mm long Shannon reamer - Helical flute L 30 mm
ANC843	Ø3.0 mm long Shannon Isham reamer - Straight flute L 30 mm

**The percutaneous reamers are only available on demand. They are supplied in a sterile single use package.



NEWCLIP TECHNICS
PA de la Lande Saint Martin - 45 rue des Garotières
44115 Haute Goulaine (France)
Tél. : +33 (0)2 28 21 37 12 - Fax : +33 (0)2 40 63 68 37
orders@newcliptechnics.com - www.newcliptechnics.com