

OSCILLATING SAW SMX 0001 SMX 0002

IDEAL FOR USE IN ORTHOPAEDIC OSTEOTOMIES AND OTHER BONE-CUTTING PROCEDURES. THE OSCILLATING SAW PROVIDES CONSTANT POWER FOR THE ENTIRE PROCEDURE.

ADDITIONAL ACCESSORIES, EXCEPT THE BLADE, ARE NOT REQUIRED. THE OSCILLATING SAW IS LIGHTWEIGHT, ALLOWING FOR OPTIMAL MANOEUVRABILITY, MINIMIZING WRIST FATIGUE. THE LINVATEC OR STRYKER TYPE LOCKING HUB OF THE SAW ACCEPTS STANDARD OSCILLATING SAW BLADES. THE OSCILLATING SAW COMES IN A STERILE PACKAGE, PROVIDING THE USER WITH THE CONFIDENCE OF STERILITY AND AVAILABILITY WHEN REQUIRED.

- 0-12,500 CYCLES PER MIN
- WEIGHT 600 GRAMS
- CUTTING ARC 4.5°

RECIPROCATING SAW SMX 0004

DESIGNED FOR USE WHEN PERFORMING VERTICAL OSTEOTOMIES SUCH AS CUTTING THE HEAD/NECK OF THE FEMUR OFF DURING TOTAL HIP PROCEDURES, UNI CONDYLAR KNEE REPLACEMENTS, GENERAL ORTHOPAEDIC OSTEOTOMIES, AND OTHER BONE CUTTING PROCEDURES. THE SOLOMAX RECIPROCATING SAW COMES PRE-STERILISED AND ACCEPTS STANDARD LINVATEC STYLE RECIPROCATING SAW BLADES.

- SPEED 0-12,500 CYCLE PER MIN
- BLADE STROKE 3.0MM
- WEIGHT 600 GRAMS

STERNUM SAW SMX 0008

THE SOLOMAX STERNUM SAW IS DESIGNED FOR USE IN CARDIO - THORACIC SURGERY. FOR EASE OF USE, THIS VERSATILE HANDPIECE IS LIGHTWEIGHT ONLY 600 GRAMS AND IS FACTORY FITTED WITH A STERNUM GUARD FOR PATIENT SAFETY AND CONVENIENCE. THE DEVICE OFFERS HIGH PERFORMANCE, RELIABILITY AND IS PRESTERILISED READY-TO-USE.

- SPEED 12500 CYCLES PER MINUTE
- BLADE STROKE 3 MM
- WEIGHT 600 GRAMS

REAMER SMX 0012 - FACTORY FITTED KEYLESS CHUCK END SMX 0014 - FACTORY FITTED ZIMMER/HUDSON END

THE SOLOMAX REAMER IS FOR USE IN ORTHOPAEDIC PROCEDURES WHICH REQUIRE DRILLING AND REAMING. THE SOLOMAX REAMER IS BEST FOR USE IN DIFFICULT PROCEDURES WHERE CONTINUOUS POWER IS CRITICAL. THE REAMER INCLUDES A 4 MM CANNULATION CHANNEL WHEN USING PINS AND GUIDE WIRES.

THE DEVICE OFFERS HIGH PERFORMANCE, RELIABILITY AND IS PRE-STERILISED READY-TO-USE. OTHER ADAPTERS ARE COMPATIBLE WITH THE SOLOMAX REAMER.

- TORQUE 6.0 NM
- ROTATION SPEED 250 RPM, FORWARD AND REVERSE
- WEIGHT 800 GRAMS


DRILL SMX 0018 - FACTORY FITTED KEYLESS CHUCK END SMX 0019 - FACTORY FITTED ZIMMER/HUDSON END

THE SOLOMAX DRILL IS DESIGNED FOR DRILLING AND SCREWING IN A VARIETY OF RECONSTRUCTIVE AND TRAUMA PROCEDURES.

FOR PINS AND GUIDE WIRES, A 4MM CANNULATED CHANNEL IS PROVIDED.

THE DEVICE IS SUPPLIED WITH A KEYLESS CHUCK, IS PRE-STERILISED, AND READY-TO-USE.

- DRILL TORQUE 1.8 NM SCREW TORQUE 1.8 NM
- DRILL SPEED 1000 RPM SCREW SPEED 250 RPM
- FORWARD & REVERSE ROTATION
- WEIGHT 700 GRAMS



POWER CONSOLE SMX 0034

THE SOLOMAX POWER CONSOLE IS AN ELECTRIC POWER SOURCE FOR THE OPERATION OF THE HANDHELD SOLOMAX POWER TOOLS. THE POWER CONSOLE FEATURES 3 CONNECTORS LOCATED ON THE FRONT PANEL, ALLOWING 3 DEVICES TO BE ATTACHED AT ONE TIME. SOLOMAX POWER TOOLS CAN BE USED EITHER WITH A POWER CONSOLE OR A RECHARGEABLE BATTERY.

- LOW VOLTAGE OUTPUT 12V/15A OUTPUT AC / WATTS / ETC
- PEAK CURRENT PROTECTION AT 25A
- ADVANCED / UNIVERSAL SWITCH MODE ARCHITECTURE
- FULLY APPROVED TO: BSEN 60601/1 / U.L / CF APPROVED
- WEIGHT 2.4 kg
- DIMENSIONS 350 MM (W) x 230 MM (D) x 90 MM (H)

BATTERY CHARGER SMX 0022

- 3 STATIONS CHARGING
- 6 HOUR CHARGING
- AUTO CUT OFF-FULL CHARGED
- AC 100V-240V VARIABLE VOLTAGE INPUT
- L.E.D DISPLAY
- WEIGHT 1.5 kg
- DIMENSIONS 380mm (W) x 90mm (H) x 220mm (D)

BATTERY PACK SMX 0026

- RECHARGEABLE CELL SET
- HIGH POWER DENSITY NI-MH CELLS
- UNIVERSAL INPUT CHARGER
- DC 12V/15A OUTPUT
- CURRENT CAPACITY AT 2.6A
- WEIGHT 650g
- DIMENSIONS 70mm (W) x 55mm (H) x 100mm (D)